

IDMS[™] INTEGRATED DATA MANAGEMENT SYSTEM WITH ILINK DATA FLOW

Data Security and Privacy White Paper





iDMS[™] INTEGRATED DATA MANAGEMENT SYSTEM WITH ILINK DATA FLOW



Introduction & Background

As medical devices evolve to incorporate more software code, dedicated operating systems and networking devices, they may also be exposed to the same cyber-threats as conventional IT systems. For example, medical devices can be infected by malware or hacked, which can impact care delivery, potentially interfere with specified control, or lead to the breach of sensitive healthcare information. With threats to the security of devices increasing, governments have enacted legislation to criminalize many of these cyberattacks and protect identifiable health information.

AdvaMed Medical Device Cybersecurity Foundational Principles provide a template for addressing cybersecurity throughout the product lifecycle, including during the design, development, production, distribution, deployment, maintenance and disposal of the device and associated data.

To address all inherent risks, manufacturers must appropriately employ a multi-modal strategy leveraging people, process and technology – both home-grown or outsourced to industry experts.

It is through this multi-modal approach that Accuray products can address strict standards for security while ensuring the expected standards in patient care and ease-of-use.

Accuray Product Security: Focus on iDMS[™]-Enabled Systems

Accuray has centralized data-management of its treatment delivery systems and software within the iDMS[™] Integrated Data Management System. This is a single, secure data-server that allows customers to easily access and share patient treatment data across any Accuray radiation therapy platform. At the same time, this centralized architecture has a robust set of security features to prevent unauthorized access while enabling patient privacy and customer ease-of-use.

As a centralized repository, iDMS-enabled products are designed with the following security measures.



SECURITY

- All computers include a white list anti-virus solution that only allows pre-defined, signed applications to operate. This protects the system against viruses and operating system vulnerabilities that could permit the deployment of viruses or malware. White list anti-virus software does not require regular updates. It also ensures that only compatible versions of Accuray software are used.
- All computer systems associated with our devices can be placed behind the included firewall. Accuray Precision[®] workstations may also (by customer request) be placed outside the firewall on the customer's network.
- Accuray uses network segmentation for performance and security purposes within the treatment system, and between the customer's network and the product network primarily for security purposes.

AUTHENTICATION/AUTHORIZATION

- Accuray systems include site-defined users and passwords. This mechanism is complemented with role-based access including a tool to define permissions per role. Patient/System data can only be accessed after the appropriate authentication has occurred. Users can view actions taken within the system within the patient and system audit logs.
- **Passwords** are maintained in the iDMS[™] database in an encrypted format at rest and in-motion. When storing encrypted passwords in iDMS, the passwords are **hashed** to protect against an attack on the database (whereby knowing one password would allow a hacker to compromise all passwords on the system).
- iDMS supports an **authentication** model that protects against repeated failed attempts to login to the system. Users that exceed the number of failed authentication attempts are temporarily disabled.

DATA INTEGRITY

- Accuray employs encryption technologies to protect sensitive data: iDMS[™] and Precision[®] encrypt all data on the file system. While Radixact[®] encrypts all patient data and CyberKnife[®] anonymizes all data on disk.
- Data is protected on the iDMS server by use of hard drives configured in a RAID 6 array. This ensures that up to 2 hard drives per array can fail prior to losing any data. Furthermore, Accuray remote monitoring can proactively warn of disk failures.
- iDMS provides file-level data integrity checking, preventing a corrupted or infected data file from being used by the system.
- Accuray provides Windows OS updates to critical systems, such as the treatment delivery system, with major software upgrades.

PRIVACY: PERSONAL DATA PROTECTION

For iDMS-enabled systems, Accuray has applied privacy-by-design principles that enable deidentification of patients' health information. Additionally, the control to transfer patient data to Accuray is at the discretion of the hospital. In the event a user intends on transferring data to Accuray, the site must manually place patient data into predefined system directories that allow Accuray personnel to access the needed information. Finally, encryption features for data at rest have been included in iDMS-enabled systems:

- Data on treatment systems is encrypted and can be de-identified in such a fashion as to remove any direct identifiers.
- When backing up patient data to a customer-provided network location, the data is automatically encrypted.
- Accuray protects PHI during service events of iDMS-enabled systems: 1) patient data if required, is de-identified prior to upload 2) log file data does not contain patient direct identifiers.

Remote Connectivity through the IoT (Internet of Things)

Accuray systems share information in real time with Accuray support teams. These include crucial systemperformance data that help prevent system failures. As well as 1-on-1 support sessions with Accuray clinical and technical support.

Remote Connectivity is managed through Axeda[®] Machine Cloud Service, a powerful, industry-leading data management technology. This secure, encrypted connection includes explicit authorization and authentication controls to ensure secure access to all information during service sessions. Axeda has attained ISO 27001:2005 certification, supporting their focus on delivering the highest levels of security, performance and availability of the Axeda Machine Cloud Service.*

iLink is an Accuray-built product that interfaces with Axeda components to establish a service connection with Accuray systems to enable remote machine diagnostics.

Data is transmitted through encrypted TomoLink/iLink conduits and kept secure using:

- 1) Transport layer security (TLS) at the communications level.
- 2) Advanced encryption standard (AES) 128 algorithm to further secure messages, and
- 3) The RSA 2048 algorithm for key data exchanges.

For real-time technical support and remote access, Accuray has partnered with LOGMEIN.

LOGMEIN's GoToAssist[®] is an on-demand application, enabling remote service sessions with Accuray support teams: Service Engineers, Clinical Applications and Physics. While their GoToMyPC[®] Corporate product is an optional remote desktop application that allows users to access Accuray treatment planning systems and data remotely within an organization's secure network.

iLink includes automatic collection of:

- System alarm events (i.e. hard drive failure, server temperature, etc.)
- Machine procedure data files including: → System sensor measurements (i.e. gun voltage, coolant flow, injector current, etc.)
 - → System workflow time measurements (i.e. procedure begin, warm-up begin, etc.)
 - → Error messages (i.e., procedure interrupted)
- System log files containing: → Software diagnostic messages → System error information
- Additionally, this secure connection can be used to:

 → Manually upload data files
 - to Accuray
 - → Remotely access systems using a site approved remote desktop sharing sessions

LOGMEIN has set an industry standard with adoption of the ISO 27001 security framework, development of tools with embedded security controls and processes that are focused on cybersecurity. Both applications employ SSL technology and 128-bit AES encryption.



Accuray iDMS[™] with iLink Data Flow

CUSTOMER RESPONSIBILITIES

- Customer is responsible for putting in place appropriate internal controls to manage system passwords and user access. This should include password strength requirements and controls for disabling access to exiting employees.
- Accuray software products should be used only on computers and networks that are properly secured in accordance with Accuray product documentation, service agreements, and instructions for use.
- Accuray systems and accessories do not permit third-party software installation by the customer (e.g., anti-virus scanners, office productivity tools, system patches, on-platform firewalls, etc.) unless documented by Accuray as an operating specification or prior written consent is attained. Unauthorized modifications to Accuray products could void warranty and alter the regulatory status of the device. Any resulting service required from unauthorized modifications is not covered under our service agreements. Such unauthorized modifications can affect the performance or safety of your device in unpredictable ways. Accuray is not responsible for equipment that has been subject to unauthorized modification.
- Accuray recommends users backup data on the hospital network so that in the case of a data loss event, a solid backup for restoration by Accuray Service personnel is available
- Accuray does not encrypt network data transmissions. It is the customer's responsibility to protect network and data transmission infrastructure from internal and external threats. This includes maintaining data security requirements and point-to-point encryption among the Accuray System firewall and the following potential data destinations or access points: The dedicated Virtual Local Area Network (VLAN), Accuray Precision[®] System workstations, Wide Area Network (WAN) Connected Treatment System, the redundant iDMS System installed on a facility Local Area Network (LAN) or WAN

Conclusion

Responsible medical device companies assess, mitigate, and constantly monitor the ever-present cybersecurity threats to critical assets. Medical device manufacturers, healthcare service providers, patients, and physicians are collectively responsible for collaborating to create a resilient, secure approach that embraces technological innovation while mitigating its associated risks. With increased information sharing, constant monitoring, and an informed understanding of the threats they face, medical device manufacturers can assess potential vulnerabilities and identify risk mitigation strategies that will ultimately strengthen security. Accuray continues to examine and evolve existing products to best accommodate the requirements of Accuray's security-minded customers while enabling a 'Patient-First' clinical practice.



Additional Information

Considering the increased focus on medical device security and compliance with the HIPAA Security Rule in the US, the Healthcare Information and Management Systems Society (HIMSS) created a standard "Manufacturer Disclosure Statement for Medical Device Security" (MDS2). The MDS2 is intended to supply healthcare providers with important information that can assist them in assessing and managing the vulnerabilities and risks associated with the electronic Protected Health Information (ePHI) that is created, transmitted, or maintained by medical devices. Accuray MDS2 forms are available to customers upon request.

Axeda Machine Cloud Service

Link to Axeda White Paper "Providing Secure Remote Service and Support": http://blog.axeda.com/hs-fs/hub/514/file-13177011-pdf/docs/axeda_securitywp.pdf

LOGMEIN Resources

HIPAA Compliance Guide https://www.gotomypc.com/downloads/pdf/m/GoToMyPC_Corporate_HIPAA_Compliance_Guide.pdf

Corporate Security White Paper https://www.gotomypc.com/downloads/pdf/m/GoToMyPC_Corporate_Security_White_Paper.pdf

GoToAssist: "Remote Support Security for the Modern Enterprise" https://assets.cdngetgo.com/18/f0/2bf0fd0b4c9b8cfff91eb6addf38/gotoassist-corporate-security-whitepaper.pdf

Platform specific IT guides for Accuray CyberKnife[®], TomoTherapy[®] and Radixact[®] Systems can be requested through support@accuray.com or by calling 1.877.668.8667

All 3rd party product and company names are trademarks M or registered * trademarks of their respective holders. Use of them does not imply any affiliation with or endorsement by them.

Notes





UNITED STATES

Accuray Corporate Headquarters

1310 Chesapeake Terrace Sunnyvale, CA 94089 USA Tel: + 1.408.716.4600 Toll Free: 1.888.522.3740 Fax: + 1.408.716.4601 Email: sales@accuray.com

ASIA

Accuray Incorporated

1240 Deming Way

Madison, WI 53717

Tel: +1.608.824.2800

Fax: +1.608.824.2996

Accuray Japan K.K. Shin Otemachi Building 7F 2-2-1 Otemachi, Chiyoda-ku Tokyo 100-0004 Japan Tel: +81.3.6265.1526

Tel: +81.3.6265.1526 Fax: +81.3.3272.6166

Accuray Asia Ltd. Units 910-11, Ocean Centre Harbour City 5 Canton Road, T.S.T Hong Kong Tel : +852.2247.8688 Fax : +852.2175.5799

Accuray Accelerator Technology (Chengdu) Co., Ltd. No. 8, Kexin Road Hi-Tech Zone (West Area) Chengdu 611731 Sichuan

EUROPE

Accuray International Sarl Route de la Longeraie 9 (3rd floor) CH - 1110 Morges Switzerland Tel: +41.21.545.9500 Fax: +41.21.545.9501

Important Safety Information:

Most side effects of radiotherapy, including radiotherapy delivered with Accuray systems, are mild and temporary, often involving fatigue, nausea, and skin irritation. Side effects can be severe, however, leading to pain, alterations in normal body functions (for example, urinary or salivary function), deterioration of quality of life, permanent injury, and even death. Side effects can occur during or shortly after radiation treatment or in the months and years following radiation. The nature and severity of side effects depend on many factors, including the size and location of the treated tumor, the treatment technique (for example, the radiation dose), and the patient's general medical condition, to name a few. For more details about the side effects of your radiation therapy, and to see if treatment with an Accuray product is right for you, ask your doctor.

© 2018 Accuray Incorporated, All Rights Reserved. The stylized Accuray logo, CyberKnife, VSI, M6, TomoTherapy, H Series, Tomo, TomoH, TomoHD, TomoHDA, TomoEDGE, TomoHelical, TomoDirect, Hi Art, PlanTouch, PreciseART, PreciseRTX, Radixact, Accuray Precision, iDMS, Iris, Xchange, RoboCouch, InCise, MultiPlan, Xsight, Synchrony, InTempo, TxView, PlanTouch, QuickPlan, TomoHelical, TomoEDGE, CTrue, VoLO, Planned Adaptive, TQA, TomoLink, TomoPortal, OIS Connect and AERO are trademarks or registered trademarks of Accuray Incorporated, in the United States and other countries and may not be used or distributed without written authorization from Accuray Incorporated. Use of Accuray Incorporated's trademarks requires written authorization from Accuray Incorporated. MKI-TxPIg-0318-0055